

2012 INFANT MANIKIN COMPARISON CHART

	<p>PP-IM-100 INFANT CPR MANIKIN</p> <p>PRESTAN</p> <ul style="list-style-type: none"> • Realistic looking • Simple lung bag insertion • Secured face shields • Visible chest rise • Head tilt & chin lift 	<p>MSRP 99.00</p> <ul style="list-style-type: none"> • Optional CPR Monitor on diaper allows for instant visual and audible feedback • Easily cleaned • FREE CPR Mask & Face Shields! • FREE Shipping at CPR Savers! 	<p>Dimensions: 10"x6"x22" Weight: 6 lbs. Warranty: 3 years 100% Latex Free</p>
	<p>LF06002U TMAN2 INFANT</p> <p>CPR PROMPT</p> <ul style="list-style-type: none"> • Race neutral in blue or realistic in tan color • Visible chest rise • Removeable chest plate for easy lung bag insertion 	<p>MSRP 59.95</p> <ul style="list-style-type: none"> • Inner chest plate gives "bottoming out" feel to ensure proper technique • Shoulder blade landmarks • Features Adam's apple & carotid arteries • Pajama markings for proper hand positions 	<p>Weight: 1.3 lbs. Warranty: 1 year 100% Latex Free</p>
	<p>LF03720U BABY BUDDY</p> <p>Nasco Life/form</p> <ul style="list-style-type: none"> • Sanitary one-piece lung, mouth protection system • Head tilt • Visible chest rise • Palpable landmarks 	<p>MSRP 90.50</p> <ul style="list-style-type: none"> • Easy-to-use, lightweight & compact • Removeable chest for easy lung/mouth bag installation • Can only be used with AED Trainer Pads made for foam manikins 	<p>Dimensions: 14"x7"x20" Weight: 3 lbs. Warranty: 3 year 100% Latex Free</p>
	<p>50000 BABY ANNE</p> <p>Laerdal helping save lives</p> <ul style="list-style-type: none"> • Extreme realism • Visible chest rise • Head tilt & chin lift 	<p>MSRP 96.03</p> <ul style="list-style-type: none"> • Oral and nasal passages allow realistic nose pinch • Natural obstruction of the airway • Foreign-body airway obstruction feature: practice back-blows and chest-thrusts 	<p>Dimensions: 8"x6"x 20" Weight: 3.3 lbs. Warranty: 1 year</p>