












# 2012 MANIKIN COMPARISON CHART

	<p>PP-AM-100 <b>ADULT CPR MANIKIN</b></p>  <ul style="list-style-type: none"> <li>• Realistic looking</li> <li>• Simple lung bag insertion</li> <li>• Secured face shields</li> <li>• Head tilt &amp; chin lift</li> <li>• Realistic landmarks</li> </ul>	<p><b>MSRP 104.75</b></p> <ul style="list-style-type: none"> <li>• Optional CPR Monitor allows for visual feedback</li> <li>• Realistic nose pinch</li> <li>• FREE CPR Mask &amp; Face Shields!</li> <li>• <b>FREE Shipping at CPR Savers!</b></li> </ul>	<p>Dimensions: 14"x7"x23" Weight: 7 lbs. Warranty: 1 year  100% Latex Free</p>
	<p>LF06001U <b>TMAN1 ADULT-CHILD</b></p>  <ul style="list-style-type: none"> <li>• Race neutral in blue or realistic in tan color</li> <li>• Visible chest rise</li> <li>• Head tilt &amp; chin lift</li> <li>• Proper landmarks</li> </ul>	<p><b>MSRP 84.95</b></p> <ul style="list-style-type: none"> <li>• Clicker allows for instant audible feedback</li> <li>• Age Selector Dial: Adjusts compression depth for Adult &amp; Child CPR options!</li> <li>• Can only be used with AED Trainer Pads made for foam manikins</li> </ul>	<p>Dimensions: 14"x7"x20" Weight: 2.4 lbs. Warranty 1 year  100% Latex Free</p>
	<p>LF03693U <b>BASIC BUDDY</b></p>  <ul style="list-style-type: none"> <li>• Sanitary airway system</li> <li>• Visible chest rise</li> <li>• Head tilt &amp; chin lift</li> <li>• Proper landmarks</li> </ul>	<p><b>MSRP 90.50</b></p> <ul style="list-style-type: none"> <li>• Lightweight &amp; compact</li> <li>• Age Selector Dial: Adjusts compression depth for Adult &amp; Child CPR options!</li> <li>• Can only be used with AED Trainer Pads made for foam manikins</li> </ul>	<p>Dimensions: 19"x16"x22" Weight: 3 lbs. Warranty: 1 year  100% Latex Free</p>
	<p>2002001 <b>LITTLE ANNE</b></p>  <ul style="list-style-type: none"> <li>• Realistic facial features &amp; movable jaw</li> <li>• Visible chest rise</li> <li>• Head tilt &amp; chin lift</li> <li>• Realistic landmarks</li> </ul>	<p><b>MSRP 218.00</b></p> <ul style="list-style-type: none"> <li>• Optional clicker allows for instant audible feedback</li> <li>• Natural obstruction of the airway</li> <li>• Realistic nose pinch</li> </ul>	<p>Dimensions: 13.5"x8.5"x 25" Weight: 8.5 lbs. Warranty: 1 year</p>

	<p>WLPMA1 <b>CPR PRACTI-MAN</b></p> 	<p><b>MSRP 139.95</b></p> <ul style="list-style-type: none"> <li>• Extremely lightweight</li> <li>• Realistic skin-like look &amp; feel</li> <li>• Head tilt &amp; chin lift</li> <li>• Realistic landmarks</li> <li>• Clicker allows for instant audible feedback</li> <li>• Age Selector Dial: Adjusts compression depth for Adult &amp; Child CPR options!</li> <li>• Requires nose pinch to complete compressions</li> </ul>	<p>Dimensions: 14"x7.5"x16" Weight: 4.6 lbs. Warranty: 2 years  100% Latex Free</p>
	<p>2800 <b>BRAD CPR MANIKIN</b></p> 	<p><b>MSRP 180.00</b></p> <ul style="list-style-type: none"> <li>• Realistic look</li> <li>• Head tilt &amp; chin lift</li> <li>• Realistic landmarks including notch</li> <li>• Easily manipulated airway for obstruction or choking situations</li> <li>• Realistic vinyl plastic over polyurethane foam</li> <li>• Longer torso for abdominal thrusts</li> </ul>	<p>Dimensions: 18"x10"x28" Weight: 15 lbs. Warranty: 1 year  100% Latex Free</p>
	<p>A259004000 <b>AMBU CPR PAL</b></p> 	<p><b>MSRP 273.60</b></p> <ul style="list-style-type: none"> <li>• Lightweight &amp; compact</li> <li>• Visible chest rise</li> <li>• Head tilt &amp; chin lift</li> <li>• Realistic landmarks</li> <li>• Stomach inflation indicated by whistling sound</li> <li>• Airways only open when head is properly hyper extended</li> </ul>	<p>Dimensions: 15"x9.8"x28" Weight: 4 lbs. Warranty: 1 year</p>
	<p>A268003000 <b>AMBU SAM</b></p> 	<p><b>MSRP 265.00</b></p> <ul style="list-style-type: none"> <li>• Sanitary airway system</li> <li>• Visible chest rise</li> <li>• Head tilt &amp; chin lift</li> <li>• Realistic landmarks</li> <li>• Feel ventilation through nose &amp; mouth</li> <li>• Store the head inside the body for easy storage &amp; transportation</li> <li>• Compression indicator confirms correct procedure</li> </ul>	<p>Dimensions: 12"x7.9"x27" Weight: 4 lbs. Warranty: 1 year  100% Latex Free</p>



LF03716U

## BASIC CPARLENE

**MSRP 595.00**

**CPARLENE®**

- Fully articulated head, neck, & jaw
- Visible chest rise
- Head tilt & chin lift
- Realistic landmarks
- Obstructed airways capability
- Sanitary airway system
- Simulated carotid pulse
- Upgradable & interchangeable accessories & add-ons, no conversion kit required

Dimensions: 13.5"x8.5"x25"  
Weight: 29 lbs.  
Warranty: 3 years  
Contains dry natural rubber



S308

## CPR SIMON

**MSRP 625.00**

**Gaumard®**

- Fully articulating head
- Extreme facial realism
- Visible chest rise
- Head tilt & chin lift
- Realistic landmarks
- Clicker allows for instant audible feedback
- Realistic inside & out: Contains ribcage, lungs, & heart
- Obstructed airway when head is forward

Weight: 7 lbs.  
Warranty: 1 year

## F.A.Q.

### WHICH MANIKIN IS BEST?

All the CPR-AED Training Manikins we carry at CPR Savers & First Aid Supply LLC will teach you the proper techniques for performing CPR in an emergency situation! There is no right or wrong manikin to choose. Some of our manikins are hyper realistic, some are extremely lightweight, and some are more advanced than others. It depends on your budget and your purpose! If you are looking to purchase manikins for an entire class, we carry packs of manikins at great prices in our store! [www.CPR-Savers.com](http://www.CPR-Savers.com)

### WHAT ARE FACE SHIELDS AND WHY DO I NEED THEM?

Face shields are single use, sanitary barriers that you may use during CPR Training. They're a great way to share manikins during class without the worry of cross-contamination!

### CAN I USE AN AED TRAINER ON MY MANIKIN?

AED Trainers can be used on most of these manikins. You must use special pads for foam manikins.

### WHAT IS THE DIFFERENCE IN INFANT, CHILD, & ADULT MANIKINS?

They require different amounts of pressure when performing compressions. Some of the manikins we carry have an Adult/Child Switch, which is a great money saver!

### ARE THESE ALL THE MANIKINS AVAILABLE?

No. These are only the popular, single, torso-only models. We have a huge selection in our online store.

### WHAT ARE CPR MONITORS?

Many of the manikins make an audible clicking sound that will confirm you are doing the compressions correctly. Some manikins have lights that will give you immediate visual feedback.

### WHAT IS MSRP?

Manufacturer Suggested Retail Price. The actual price may be different than this.