



Body Loc

Product#: 30030

Includes:

- Carry Case: 14" L x 6" W x 8" H
- 1 Black/Orange Pelvis Strap (78" L x 3" W)
- 1 Black/Yellow Femur (73" L x 3" W)
- 1 Black/Red Chest Strap (68" L x 3" W)
- 1 Black/Green Upper Tibia Strap (63" L x 3" W)
- 1 Black/Blue Lower Tibia Strap (58" L x 3" W)

Our Body-Loc strapping system is perfect for first responders in the public and industrial markets. This original design provides extremely efficient immobilization using color coded straps with reflective anatomical positioning labels for fast, foolproof assembly. Hook and loop closures allow quick, secure connections while allowing adjustment to accommodate patient of any size.



*Kyoto & Montreal Compliant.

© Copyright 2013. Iron Duck, Built for Life, and Iron Duck logo are registered trademarks of Fleming Industries Inc. All Rights Reserved.

BUILT IN AMERICA