

IDEA | Extrication Apparatus

Product#: 35960

The Iron Duck Extrication Apparatus (IDEA) is design to help extricate patients trapped in confined areas while reducing the risk of further injury. You can use our IDEA for both spine and thorax immobilization or as a hip and pelvic splint.

The body straps are color coded for you to be able to instantly identify where each strap goes in a hurry. The IDEA is heat sealed and made with universal precautions material allowing for complete disinfection. When using our IDEA please follow all standard extrication device procedures.



*Kyoto & Montreal Compliant.

© Copyright 2013. Iron Duck, Built for Life, and Iron Duck logo are registered trademarks of Fleming Industries Inc. All Rights Reserved.

BUILT IN AMERICA