

## Oxygen Bag 2-Pockets 36002D-2PK

> Dimensions: 21"L x 14"W x 8.5"H ~ 2,499 cu in

This clamshell designed oxygen bag is designed for fast deployment of a D or Jumbo D tank, regulator and accessories. There are two full length exterior pockets to provide added storage. Zip open the front full length pocket to find a row of elastic tacking for tools and accessories of multiple sizes.

The rear full length pocket is open allowing for storage of bulky or larger size items you may need. Pull open the heavy duty zippers the full length of the bag to access the main compartment. The removable liner is made from high impact foam to help protect your tank and regulator.

A padded regulator pocket and adjustable tie-downs hold the tank securely in place during your busy day. On each side of the main compartment you will find several pockets and elastic tool holders for masks, tubes and accessories. The bottom of this oxygen bag is lined with an extra layer of UP fabric for extra durability and easy cleaning. Carry your oxygen bag with heavy duty wrap handles or with our ergonomically designed shoulder strap.

